Amazon Drop Test Evaluation for Flexible Packages

Date:

Vendor/Seller Name: ASIN:

Item: Package Weight (lbs) cannot exceed 7 lbs:

Certification (choose only one): FFP* SIOC Package Dimension (LxWxHinches):

*FFP requirements: 100% curbside recyclable packaging , Easy-To-Open, all contents removable within 120 seconds

Sealing Method: Heat Seal Adhesive Tape Resealable Tamper Evident Seal (required if resealable)

Slack is less than 2 inches

Exempt from <2inch slack (rigid and padded mailers)

No protrusions

Hazard Drop Box* used

*Hazard Drop Box must be 9 x 6 x 3 inches and 3 lbs (or equivalent)

Material Thickness \geq 2 mils LDPE or minimum tear strength of \geq 130 gram Seal strength of \geq 5.7lb/in (1 N/mm) (Adhesive Sealing Only) No holes larger than 0.25" (6.4mm)


No cosmetic or functional damage to product after testing

TEST PICTURES

	TEST TICTORES	
Pre-test 1-2-5	Pre-test 3-4-6	Post-test - BOM
110 (63) 12 3	F16-test 3-4-0	Post-test - Bolyi
Doct test. Photo of item	Doct tost. Contants marked as Posselable	Photo of Honord Dress Roywood
Post-test - Photo of item	Post-test - Contents marked as Recyclable	Photo of Hazard Drop Box used

Drop Number	Height	Orientation
1	18 in	Edge 3-4
2	18 in	Edge 3-6
3	18 in	Edge 4-6
4	18 in	Corner 3-4-6
5	18 in	Corner 2-3-5
6	18 in	Edge 2-3
7	18 in	Edge 1-2
8	36 in	Face 3
9	18 in	Face 3
10	18 in	Edge 3-4
11	18 in	Edge 3-6
12	18 in	Edge 1-5
13	18 in	Corner 3-4-6
14	18 in	Corner 1-2-6
15	18 in	Corner 1-4-5
16	36 in	Face 6
17	18 in	Face 3

Hazard Box Test: Drop the Hazard box from a height of 36 inches, and orient it such that the corner impacts the test sample's bag. Repeat the drop twice more. Total number of drops is three.



A test from an ISTA certified lab is required if it is categorized as any of the following:

- Fragile
- TV/Monitor
- · Contains liquid

For more information please visit: https://www.aboutamazon.com/packaging